



# Bead Ruptor 4

MINI BEAD MILL HOMOGENIZER

***Powerful performance  
in a compact design***

*Built upon the innovative  
design and unsurpassed  
quality of the Omni  
Bead Ruptor family.*

# Bead Ruptor 4

## MINI BEAD MILL HOMOGENIZER

ORDER NO. 25-010



### Processing Range:

25 µL to 6 mL

### Power Rating:

100-230 VAC/50-60 Hz

### Dimensions:

8.5" W x 11.5" D x 10" H

**Weight:** 16.5 lbs. (7.48 kg)

**Run Time:** 1 - 5 mins.

**Programs:** 1

**Cycles:** 1

### Performance Range:

30 x g - 150 x g

### Efficiently lyse:

- > Human tissues
- > Animal tissues
- > Plant & food
- > Micro-organisms

### Designed for personal use

- > Broad performance range — gentle mixing to high force homogenization
- > Small footprint — requires minimal bench space
- > Holds 4 x 0.5 mL / 1.5 mL / 2 mL or 1 x 7 mL tubes

### Patent Pending Snap In Tube Carriage



### Omni Bead Ruptor 4 is ideally suited for:

- > Homogenizing up to 4 samples in seconds
- > Reproducible extraction of proteins and nucleic acids
- > Simultaneously processing 4 x 0.5 mL, 4 x 1.5 mL, 4 x 2 mL or 1 x 7 mL samples



*For best results, use with Omni pre-filled bead kits*

Bio-Budget Technologies GmbH  
Dießemer Bruch 150, 47805 Krefeld  
TEL 02151 - 65 20 830

**Bio-Budget is a partner of Biolabproducts GmbH**



> **For more information —**  
**[www.bio-budget.com](http://www.bio-budget.com)**

**Bio-Budget**  
Technologies GmbH

